


BITPLANE

an **Oxford Instruments** company



“XT Cheat Sheet”

```
% get the application object
if isa(almarisApplicationID, 'Imaris.IApplicationPrxHelper')
    % called from workspace
    vlmarisApplication = almarisApplicationID;
else
    % connect to Imaris interface
    javaaddpath ImarisLib.jar
    vlmarisLib = ImarisLib;
    if ischar(almarisApplicationID)
        almarisApplicationID = round(str2double(almarisApplicationID));
    end
    vlmarisApplication = vlmarisLib.GetApplication(almarisApplicationID);
end
```

To launch XTension for debugging in Matlab type “NAMEofXTension 0”

- ‘0’ represents the Imaris application ID (by default this starts at 0)
 - if debugging fails to recognize Imaris ID, the ImarisID is likely not zero
 - close Imaris, let rest for 10-15sec and restart to establish ID=0
- no spaces in .m file are allowed

Listing Spots Objects and Selecting One

```
vSpots = vImarisApplication.GetFactory.ToSpots(vImarisApplication.GetSurpassSelection);
%SelectSpotsObject
vNumberOfSpots = 0;
vSpotsList{vSurpassScene.GetNumberOfChildren} = [];
vNamesList{vSurpassScene.GetNumberOfChildren} = [];
for vChildIndex = 1:vSurpassScene.GetNumberOfChildren
    vDataItem = vSurpassScene.GetChild(vChildIndex - 1);
    if vImarisApplication.GetFactory.IsSpots(vDataItem)
        vNumberOfSpots = vNumberOfSpots+1;
        vSpotsList{vNumberOfSpots} = vImarisApplication.GetFactory.ToSpots(vDataItem);
        vNamesList{vNumberOfSpots} = char(vDataItem.GetName);
    end
end
%test if there are any Spots Objects
if vNumberOfSpots==0
    msgbox('Please create at a spots object!');
    return;
end
vNamesList = vNamesList(1:vNumberOfSpots);
%input dialog to list all Spot objects and select one
if vNumberOfSpots > 1
    [vPair, vOk] = listdlg('ListString',vNamesList,'SelectionMode','single',...
        'ListSize',[250 150],'Name','Spots','InitialValue',1,...
        'PromptString',{'Please select spots object:'});
    if vOk<1, return, end
    vSpots = vSpotsList{vPair(1)};
else
    vSpots = vSpotsList{1};
end
```


Position Information in Volume Conversion of μm to Pixel Position

13	14	15	16	17	18
7	8	9	10	11	12
1	2	3	4	5	6

Ymax	13	14	15	16	17	18
	7	8	9	10	11	12
Ymin	1	2	3	4	5	6
	Xmin					Xmax

- numbers indicate the index in a 1DSlice in 2D
- volume size in 2D
 - 6x3 (XY)

```
%Create 2D map of XY positions in volume space
%a way of identifying the pixels masked in surface relative to XYZ spot or
vSliceXY=[repmat((linspace(vDataMin(1), vDataMax(1),
vDataSize(1)))',vDataSize(2),1)...
(repelem(linspace(vDataMin(2), vDataMax(2),
vDataSize(1)),vDataSize(2)))'];
%Identify the Z-position from the Slice#
vSliceZ=(linspace(vDataMin(3), vDataMax(3), vDataSize(3)))';
```

```
%This code find the closest pixel position to the sub voxel spot positions
for vIndexZ = 1:vDataSize(3)
    %Select spots that lie closest to current slice
    idx=subvoxelSpotPosition(:,3) > Step(vIndexZ) & ...
        subvoxelSpotPosition(:,3) < Step(vIndexZ+1);
    if idx==true
        CurrentSliceSpots=subvoxelSpotPosition(idx,:);
        %Find the pixel row and column in XY that most closely fits Spot
        position
        CountY=round((CurrentSliceSpots(:,2)-
vDataMin(2))/Yvoxelspacing,0);
        CountX=round((CurrentSliceSpots(:,1)-
vDataMin(1))/Xvoxelspacing,0);
        SpotIndexfor1D=(CountY)*vDataSize(1)+CountX;
        %Calculate the total index if assume entire Volume of slices
        SpotIndexWholeVoulme=SpotIndex+(vIndexZ*SliceInterval);
        vSlice(SpotIndexSlice1D)=255;
    end
end
```

Setting New Spots Object in Imaris Scene

```
vSpotsA = vImarisApplication.GetFactory.CreateSpots;  
vSurpassScene = vImarisApplication.GetSurpassScene;  
vSpotsA.Set(vSpotsPosXYZ, vSpotsPosT, vSpotsRadius);  
vSpotsA.SetName(sprintf('NEW SURFACE'));  
vSpotsB.SetColorRGBA(255*256*256);  
%Add Spots to scene  
vImarisApplication.GetSurpassScene.AddChild(vSpotsA, -1);
```

NOTE: this is the format for the variables (equal # rows required for each one)

```
PositionsXYZ = [10 37 10; 15 20 8];  
aIndices = [0 0]; % Indices start at zero for first time point  
aRadii = [0.5 0.5];
```

Identify All Track IDs for All Spots Tracked

```
vSpots = vImarisApplication.GetSurpassSelection;

edges = vSpots.GetTrackEdges + 1;
edges_forspots = 1:size(vSpots.GetPositionsXYZ, 1);%Get the number of spots or surfaces
edges_forspots(:) = size(edges, 1) + 1; % initialize array to fictive edge
edges_forspots(edges(:, 1)) = 1:size(edges, 1);
edges_forspots(edges(:, 2)) = 1:size(edges, 1);
trackid_foredges = [vSpots.GetTrackIds; 0]; % add fictive track id
trackid_forspots = trackid_foredges(edges_forspots);
tempx = double(trackid_forspots);
vtrackID = tempx-1000000000;%track IDs as a single integer
vtrackIDmax = max(vtrackID);%Identifies the total number of tracks in the dataset
%Loop each trackID
for trackloop = 0:vtrackIDmax;
    vSpotsT = vtrackID == trackloop;
    tracklength = sum(vSpotsT) ;
    vSpotsIndex = find(vSpotsT);%
    %Loop through each object in the TrackID
    for c = 1:tracklength-1;
        vWorkingSpotsIds=vIds(vSpotsIndex,:);
        vWorkingId=double(vWorkingSpotsIds(c));
        vWorkingTimeIndex=vSpots.GetTimeIndex(vWorkingId);
    end
end
end
```

Identify All Track IDs for All Surfaces Tracked

```
vSurfaces = vlmarisApplication.GetSurpassSelection;

edges = vSurfaces.GetTrackEdges + 1;
edges_forsurfaces = 1:size(Spots.GetPositionsXYZ, 1);%Get the number of spots or surfaces
edges_forsurfaces( : ) = size(edges, 1) + 1; % initialize array to fictive edge
edges_ forsurfaces (edges(:, 1)) = 1:size(edges, 1);
edges_ forsurfaces (edges(:, 2)) = 1:size(edges, 1);
trackid_foredges = [vSurfaces.GetTrackIds; 0]; % add fictive track id
trackid_ forsurfaces = trackid_foredges(edges_ forsurfaces);
temp_x = double(trackid_ forsurfaces);
vtrackID = temp_x-10000000000;%Track IDs as a single integer.
vtrackIDmax = max(vtrackID);%Identifies the total number of tracks in the dataset
%Loop each trackID
for trackloop = 0:vtrackIDmax;
    vSurfacesT = vtrackID == trackloop;
    tracklength = sum(vSurfacesT) ;
    vSurfacesIndex = find(vSurfacesT);%
    %Loop through each object in the TrackID do whatever
    for c = 1:tracklength;
        vWorkingSurfacesIds=vIds(vSurfacesIndex,:);%Get all surface IDs in track
        vWorkingId=double(vWorkingSurfacesIds(c));%Get first object ID in Track
        vWorkingTimeIndex=vSurfaces.GetTimeIndex(vWorkingId);

    end
end
```

```
vFilamentsIndexT = vFilaments.GetTimeIndex(FilamentIndex);  
vFilamentsXYZ = vFilaments.GetPositionsXYZ(FilamentIndex);  
vFilamentsRadius = vFilaments.GetRadii(FilamentIndex);  
vSegmentIds=unique(vFilamentsEdgesSegmentId);%Identify unique filament  
segmentID  
for vBranchIndex=1:(vNumberOfDendriteBranches)  
    %Set the ID for dendrite segment  
    wSegmentIndex = vSegmentIds(vBranchIndex);  
    %Logical argument to identify spots in segment  
    vSpotsT = vFilamentsEdgesSegmentId == wSegmentIndex;  
    vSpotsIndex = find(vSpotsT');  
    %Identify position for dendrite segment  
    %Test with new method of filtering  
    vDendriteEdgesWorking=vFilamentsEdges(vSpotsIndex,:);  
    vEdgesUnique=unique(vDendriteEdgesWorking);  
    vDendritePositionsWorking=vFilamentsXYZ(vEdgesUnique,1:3);  
    vDendriteRadiusWorking=vFilamentsRadius(vEdgesUnique,1);  
    vTypesWorking=vTypes(vEdgesUnique,:);  
end
```


Reorder Filament Points per Dendrite

%find single indices, identifying first and last points

```
vDendriteStartEndIndex=find(accumarray(vDendriteEdgesWorking(:,1))==1);
for k=1:numel(vDendritePositionsWorking(:,1))
    if k==1 %for first point
        vDendritePositionsNEW(k,:)=vFilamentsXYZ(vDendriteStartEndIndex(1),:);
        test=vDendriteEdgesWorking(:,1)==vDendriteStartEndIndex(1);
        if test==false;
            test=vDendriteEdgesWorking(:,2)==vDendriteStartEndIndex(1);
        end
        %Find next Filament Index
        vDendriteNextFilamentIndex=vDendriteEdgesWorking(test,:);
        vDendriteNextFilamentIndex(vDendriteNextFilamentIndex==vDendriteStartEndIndex(1))=[];
        %remove current edge row
        vDendriteEdgesWorking(test,:)=[];
        test(test,:)=[];

    else % each additional point
        vDendritePositionsNEW(k,:)=vFilamentsXYZ(vDendriteNextFilamentIndex,:);
        test=vDendriteEdgesWorking(:,1)==vDendriteNextFilamentIndex;
        if test==false;
            test=vDendriteEdgesWorking(:,2)==vDendriteNextFilamentIndex;
        end
        vDendriteLastFilamentIndex=vDendriteNextFilamentIndex;
        %Find next Filament Index
        vDendriteNextFilamentIndex=vDendriteEdgesWorking(test,:);
        vDendriteNextFilamentIndex(vDendriteNextFilamentIndex==vDendriteLastFilamentIndex)=[];
        %remove current edge row
        vDendriteEdgesWorking(test,:)=[];
        test(test,:)=[];
    end
end
end
```

%Where the variable NewStat represents the values to be used in the stat. They should be equal to number of surfaces, spots or tracks that are a part to the object.

%Variation #1

%adding new stat to spot/surface object for each time point in objects %already tracked in Imaris

vInd = 1:numel(NewStat);%Quantifies the number of values.

vIds = vObject.GetIds;%Gets all Ids to correspond to Imaris 8.2 Ids

vUnits(vInd) = { char(vImarisApplication.GetDataSet.GetUnit) };%gets default units (um) and converts to character

vFactors(vInd) = {'Spot'};%Generates

vFactors(2, vInd) = num2cell(vNewSpots.GetIndicesT);

vFactors(2, vInd) = cellfun(@num2str, vFactors(2, vInd), 'UniformOutput', false);

vFactorNames = {'Category', 'Time'};

vNames(vInd) = {'NewStatName'};%Name the statistic as it will appear in Imaris

NewStat=NewStat';%Transposes the values from vertical column to a horizontal arrangement

vNewObject.AddStatistics(vNames, NewStat, vUnits, vFactors, vFactorNames, vIds);

Adding a New Track Statistic

%Set New Track Statistic

vIndT=1:vtrackIDmax+1;%Total number of tracks

vldsT=0:vtrackIDmax;%Total number from tracks starting at 0

vldsT=vldsT+1000000000;%Conversion to Tracks reported by Imaris

vUnitsT(vIndT) = {'seconds'};

vFactorsT(vIndT) = {'Track'};

vNamesT(vIndT) = {' NewStat Name'};%Statistic name

vFactorNamesT = {'Category'};

vObject.AddStatistics(vNamesT, vContactTime, vUnitsT, vFactorsT, vFactorNamesT,
vldsT);

NOTE: vtrackIDmax, calculated from previous script to identify trackIDs for all
objects

Adding a New Overall Statistic

```
vInd=1:aSizeT;  
vlds(vInd)=0;%For overall stats all lds are equal to 0  
vUnits(vInd) = {'UnitName'};%{  
char(vImarisApplication.GetDataSet.GetUnit) };  
Indices=1:aSizeT;%These range for each time point  
starting at 1  
vFactors(vInd) = {'Overall'};  
vFactors(2, vInd) = num2cell(Indices);  
vFactors(2, vInd) = cellfun(@num2str, vFactors(2,  
vInd), 'UniformOutput', false);  
vFactorNames = {'Overall','Time'};  
vNames(vInd) = {sprintf('NEW stat NAME')};  
vSpots.AddStatistics(vNames,  
PercentageContactsperTimpoint', vUnits, vFactors,  
vFactorNames, vlds);
```

%Get All Statistics

```
vAllStatistics = vSpots.GetStatistics;  
vNames        = cell(vAllStatistics.mNames);  
vValues        = vAllStatistics.mValues;  
vUnits         = cell(vAllStatistics.mUnits); % not used  
vFactors       = cell(vAllStatistics.mFactors);  
vFactorNames   = cellstr(char(vAllStatistics.mFactorNames));  
vIds           = vAllStatistics.mIds;  
vObjectIndex = strmatch('Intensity Min', vNames);  
vObjectValues = vValues(vObjectIndex,:);
```



```
vRGBA=[255 0 0 0];%red
```

```
vRGBA=[0 255 0 0];%green
```

```
vRGBA=[0 0 255 0];%blue
```

```
vRGBA = uint32(vRGBA * [1; 256; 256*256; 256*256*256]);
```

```
vNewObject.SetColorRGBA(vRGBA);
```

Find Imaris Version, Quit Imaris, Start New Instance with Next Imaris ApplicationID

%Get Imaris version

```
aVersion = char(vlmarisApplication.GetVersion);  
almarisFolderEnd = strfind(aVersion, '[');  
if numel(almarisFolderEnd) ~= 1  
    msgbox('Invalid Imaris version')  
    return  
end  
almarisFolder = aVersion(1:almarisFolderEnd);  
aDelimiters = strfind(almarisFolder, '-');  
if numel(aDelimiters) == 2  
    almarisFolder(aDelimiters(2)) = [];  
    almarisFolder(aDelimiters(1)) = '';  
end
```

%Quit Imaris Application

```
vlmarisApplication.SetVisible(~vlmarisApplication.GetVisible);  
vlmarisApplication.Quit;
```

%Start new Imaris instance with ID99

```
eval(['! C:\\Program Files\\Bitplane\\', almarisFolder, '\\Imaris.exe &'])
```

%%

```
almarisApplicationID=almarisApplicationID+1;  
vlmarisApplication = vlmarisLib.GetApplication(almarisApplicationID)
```



BITPLANE

an **Oxford Instruments** company

Find us on  