




BITPLANE

an **Oxford Instruments** company

The Anatomy of an Imaris Matlab XTension



Introduction to Imaris Interface

Connect to Imaris in Java

To connect to an **Imaris** instance, it is first necessary to connect to the `ImarisServer::IServer` at port 4029 using **ICE**. The package `Ice` is included in `Ice.jar`. The package `ImarisServer` is generated compiling `bplmarisServerIce.ice`.

Example in Java:

```
ImarisServer.IServerPrx GetServer() {  
    Ice.Communicator vCommunicator = Ice.Util.initialize();  
    Ice.ObjectPrx vObject = mCommunicator.stringToProxy("ImarisServer:default -p 4029");  
    ImarisServer.IServerPrx vServer = ImarisServer.IServerPrxHelper.checkedCast(vObject);  
    return vServer;  
}
```


Connect to Imaris in Matlab

The code shown to connect to **Imaris** in Java is provided as a utility in `ImarisLib.jar` (this file is located in `imaris_installation_folder/XT/matlab`). This jar file can be used to connect to **Imaris** from Matlab, that does not directly supports **ICE**.


Example in Matlab:

```
function aImarisApplication = GetImaris  
javaaddpath ImarisLib.jar;  
vImarisLib = ImarisLib;  
vObjectId = 0; % this might be replaced by "vObjectId = <a name=getobjectId><b>GetObjectId</b></a>" (see later)  
aImarisApplication = vImarisLib.GetApplication(vObjectId);
```

As the first **Imaris** instance that registers itself to the Server is assigned to ID zero, an `aObjectId` equal to zero will work in most of the cases. `ImarisLib.jar` grants access to the Server; this can be useful in case of multiple instances of **Imaris** are started.



Simple XTension Workflow



Inserts XT option in
the 'surfaces function'
folder in the Image
Processing menu tab

Inserts XT option in
the Surpass scene
(Gear tab)

```
7  %<CustomTools>
8  %
9  % <Menu>
10 %
11 % <Submenu name="Surfaces Functions">
12 %
13 % <Item name="BiofilmAnalysis" icon="Matlab">
14 % <Command>MatlabXT::XT_MJG_BiofilmAnalysis_Complete(%i)</Command>
15 % </Item>
16 % </Submenu>
17 %
18 % </Menu>
19 %
20 % <SurpassTab>
21 % <SurpassComponent name="bpSurfaces">
22 % <Item name="BiofilmAnalysis" icon="Matlab">
23 % <Command>MatlabXT::XT_MJG_BiofilmAnalysis_Complete(%i)</Command>
24 % </Item>
25 % </SurpassComponent>
26 % </SurpassTab>
27 %
28 % </CustomTools>
```

- Inserting for other Surpass objects (GearTab)
 - bpSpots
 - bpFilaments
 - bpCells
- One XTension can be added to more than one Surpass Object (See Distance Transform)
- These are always commented out of code with “%” in Matlab

- Initiating the connection between Matlab and a specific Imaris instance

```
function XTSpotsCloseToSurface(aImarisApplicationID)


% connect to Imaris interface
if ~isa(aImarisApplicationID, 'Imaris.IApplicationPrxHelper')
    javaaddpath ImarisLib.jar
    vImarisLib = ImarisLib;
    if ischar(aImarisApplicationID)
        aImarisApplicationID = round(str2double(aImarisApplicationID));
    end
    vImarisApplication = vImarisLib.GetApplication(aImarisApplicationID);
else
    vImarisApplication = aImarisApplicationID;
end
```

- Each instance of Imaris application that is running has a unique ID (starting at 0)
- This code sets the current ID when you run XTension from Imaris.

NOTE: when debugging Imaris in Matlab, you will have to specify the ID to launch code manually

- When a single Surpass scene object is selected in Imaris (Spots, Surface, Filament or Cell)

```
61 % get the spots and the surface object
62 - vSpots = vImarisApplication.GetFactory.ToSpots(vImarisApplication.GetSurpassSelection);
63 - vSurfaces = vImarisApplication.GetFactory.ToSurfaces(vImarisApplication.GetSurpassSelection);
64
65 - vSpotsSelected = ~isequal(vSpots, []);
66 - vSurfaceSelected = ~isequal(vSurfaces, []);
67 - if vSpotsSelected
68 -     vParent = vSpots.GetParent;
69 - elseif vSurfaceSelected
70 -     vParent = vSurfaces.GetParent;
71 - else
72 -     vParent = vSurpassScene;
73 - end
```





- This code checks if any Surpass object is selected
 - i.e. attached screenshot show Spots2, as selected*
- Enters this Object as an active variable in the Matlab environment

Selecting Surpass Object(s) (Method #2)

```
69 % the user has to create a scene with some spots
70 - vSurpassScene = vImarisApplication.GetSurpassScene;
71
72 - vNumberOfSpots = 0;
73 - vSpotsList{vScene.GetNumberOfChildren} = [];
74 - vNamesList{vScene.GetNumberOfChildren} = [];
75 - for vChildIndex = 1:vScene.GetNumberOfChildren
76 -     vDataItem = vScene.GetChild(vChildIndex - 1);
77 -     if vImarisApplication.GetFactory.IsSpots(vDataItem)
78 -         vNumberOfSpots = vNumberOfSpots+1;
79 -         vSpotsList{vNumberOfSpots} = vImarisApplication.GetFactory.ToSpots(vDataItem);
80 -         vNamesList{vNumberOfSpots} = char(vDataItem.GetName);
81 -     end
82 - end
83
84 - if vNumberOfSpots<2
85 -     msgbox('Please create at least 2 spots objects!');
86 -     return;
87 - end
88 - vNamesList = vNamesList(1:vNumberOfSpots);
89 - vPair = [];
90 - while length(vPair) ~= 2
91 -     [vPair, vOk] = listdlg('ListString',vNamesList,'SelectionMode','multiple',...
92         'ListSize',[250 150],'Name','Colocalize spots','InitialValue',[1,2], ...
93         'PromptString',{'Please select the 2 spots to colocalize:'});
94 -     if vOk<1, return, end
95 -     if length(vPair) ~= 2
96 -         vHandle = msgbox(['Please select two (2) objects. Use "Control" and left ', ...
97             'click to select/unselect an object of the list.']);
98 -         uiwait(vHandle);
99 -     end
100 - end
101
102 - vSpots1 = vSpotsList{vPair(1)};
103 - vSpots2 = vSpotsList{vPair(2)};
```

 Collate and
list all Spots
Objects

 Generate Input
dialog to display
list and select
one or more
Objects

 Set selected
Objects as
active variables

- Retrieving some basic image Properties from Volume

```
2  %Get Image Data parameters
3  vDataMin = [vImarisApplication.GetDataSet.GetExtendMinX,...
4             vImarisApplication.GetDataSet.GetExtendMinY,...
5             vImarisApplication.GetDataSet.GetExtendMinZ];
6  vDataMax = [vImarisApplication.GetDataSet.GetExtendMaxX,...
7             vImarisApplication.GetDataSet.GetExtendMaxY,...
8             vImarisApplication.GetDataSet.GetExtendMaxZ];
9  vDataSize = [vImarisApplication.GetDataSet.GetSizeX,...
10             vImarisApplication.GetDataSet.GetSizeY,...
11             vImarisApplication.GetDataSet.GetSizeZ];
12  vNumberOfChannels = vImarisApplication.GetDataSet.GetSizeC;
13  Xvoxelspacing = (vDataMax(1)-vDataMin(1))/vDataSize(1);
14  Yvoxelspacing = (vDataMax(2)-vDataMin(2))/vDataSize(2);
15  Zvoxelspacing = (vDataMax(3)-vDataMin(3))/vDataSize(3);
16
17  %Get Spots information
18  vSpotsXYZ = vSpots.GetPositionsXYZ;
19  vSpotsTime = vSpots.GetIndicesT;
20  vSpotsRadius = vSpots.GetRadiiXYZ;
21  vSpotsEdges = vSpots.GetTrackEdges;
```

Getting the basic
Spot parameters
from Imaris



- Retrieving some basic Surpass object values

Surfaces

GetCenterOfMass (int aSurfaceIndex)
GetIds ()
GetMask (float aMinX, float aMinY, float aMinZ)
GetNumberOfSurfaces ()
GetSelectedIds ()
GetSelectedIndices ()
GetSingleMask (int aSurfaceIndex, float aMinX, float aMinY, float aMinZ)
GetSurfaceData (int aSurfaceIndex)
GetSurfaceDataLayout (int aSurfaceIndex)
GetSurfaceNormals (int aSurfaceIndex)
GetTimeIndex (int aSurfaceIndex)
GetTimePoint (int aIndexT)
GetTrackEdges ()
GetTrackIds ()
GetTransform (int aSurfaceIndex)

Spots

Get ()
GetIds ()
GetIndicesT ()
GetPositionsXYZ ()
GetRadii ()
GetRadiiXYZ ()
GetSelectedIds ()
GetSelectedIndices ()
GetTimePoint (int aIndexT)
GetTrackEdges ()
GetTrackIds ()

Filaments

GetBeginningVertexIndex (int aFilamentIndex)
GetColorSpinesRGBA ()
GetEdges (int aFilamentIndex)
GetEdgesSegmentId (int aFilamentIndex)
GetFilamentsList (tInts aFilamentIndices)
GetFilamentTrackEdges ()
GetFilamentTrackIds ()
GetIds ()
GetNumberOfFilaments ()
GetPositionsXYZ (int aFilamentIndex)
GetRadii (int aFilamentIndex)
GetSelectedIds ()
GetSelectedPositionIndex (int aFilamentIndex)
GetSelectedPositionIndices (int aFilamentIndex)
GetTimeIndex (int aFilamentIndex)
GetTypes (int aFilamentIndex)
GetVertexTrackEdges ()
GetVertexTrackIds ()

- Identify TrackID's and loop each track

```
33 %Identify all track IDs for all spots tracked
34 %Here is a beautiful code that does the mapping without the need of a loop
35 %(the mapping is somehow complicated because one spot might be part of
36 %different edges, usually two).
37 edges = Spots.GetTrackEdges + 1;
38 edges_forspots = 1:size(Spots.GetPositionsXYZ, 1);
39 edges_forspots( : ) = size(edges, 1) + 1; % initialize array to fictive edge
40 edges_forspots(edges(:, 1)) = 1:size(edges, 1);
41 edges_forspots(edges(:, 2)) = 1:size(edges, 1);
42 trackid_foredges = [Spots.GetTrackIds; 0]; % add fictive track id
43 trackid_forspots = trackid_foredges(edges_forspots);
44 tempx = double(trackid_forspots);
45 vtrackID = tempx-1000000000;
46 vtrackIDmax = max(vtrackID);
47
48 %Start loop for each track
49 for trackloop = 0:vtrackIDmax;
50     %Identify trackID and create working position matrix array for that trackID
51     vSpotsT = vtrackID == trackloop;
52     tracklength = sum(vSpotsT) ;
53     vSpotsIndex = find(vSpotsT');
54     vWorkingPositionXYZ = vPositionXYZ(vSpotsIndex,:);
55     vWorkingSpotsT = vSpotsT(vSpotsIndex,:);
56     vWorkingSpotsTime = vSpotsTime(vSpotsIndex,:);
57
58     % Do track processing on per track basis
59
60
61 end
```

- Loop all points for each dendrite

```
3  for FilamentIndex=0:NumberOfFilaments
4      vFilamentsXYZ = vFilaments.GetPositionsXYZ(FilamentIndex);
5      vFilamentsEdges = vFilaments.GetEdges(FilamentIndex) + 1;
6      vFilamentsEdgesSegmentId = vFilaments.GetEdgesSegmentId(FilamentIndex);
7      vSegmentIds=unique(vFilamentsEdgesSegmentId);%Identify unique filament segmentIDs
8      for vBranchIndex=1:vNumberOfDendriteBranches
9          %Set the ID for dendrite segment
10         wSegmentIndex = vSegmentIds(vBranchIndex);
11         vDendriteDistance(vBranchIndex)=0;
12         %Logical argument to identify spots in segment
13         vSpotsT = vFilamentsEdgesSegmentId == wSegmentIndex;
14         vSpotsIndex = find(vSpotsT');
15         %Identify position for dendrite segment
16         %Test with new method of filtering
17         vDendriteEdgesWorking=vFilamentsEdges(vSpotsIndex,:);
18         vEdgesUnique=unique(vDendriteEdgesWorking);
19         vDendritePositionsWorking=vFilamentsXYZ(vEdgesUnique,1:3);
20         vDendriteRadiusWorking=vFilamentsRadius(vEdgesUnique,1);
21         vTypesWorking=vTypes(vEdgesUnique,:);
22
23         %
24         % Do processing here for one dendrite at a time
25         %
26
27     end %loop for each dendrite segment
28 end%loop fo each Filament
29
```





BITPLANE

an **Oxford Instruments** company



Live Walkthrough of an XTension

- Copy Fiji/ImageJ to a folder where all users have write access
 - Optimal location – directly on the **C: drive**
 - **ImageJ/Fiji has no installer, easy to move this program, copy and paste the directory**
 - I am currently using FIJI/ImageJ 1.51p
 - Configure the path to the Fiji/ImageJ executable
 - **Win:** C:/Fiji.app/ImageJ-win64.exe
 - **Mac:** /Applications/Fiji.app/Contents/MacOS/fiji-macosx
- The “bridge” plugin will be installed automatically when you launch the first plugin from Imaris.



ImarisXT: ImageJ/FIJI link

- Embed possible Fiji/ImageJ Plugins
- Menu available, if Fiji/ImageJ configured



file input	image input	other input	
no (e.g. run macro on files only)	yes (e.g. save image to file)	no (e.g. store PCA result as file)	file output
yes (e.g. load image from file)	yes (e.g. gaussian image filter)	no (*)	image output
no	(yes) output not fed back to Imaris	no (e.g. do PCA based on some statistics)	other output

- Simple text file (see example below)

Sets location in Imaris menu


Sends volume from Imaris to Fiji

Launches plugin by specific name in Fiji

```
// <CustomTools>
//   <Menu name="Fiji">
//     <Submenu name="Process">
//       <Item name="UnsharpMask" icon="Fiji">
//         <Command>ImageJ::EEE_Process_UnsharpMask</Command>
//       </Item>
//     </Submenu>
//   </Menu>
// </CustomTools>
```


```
call("Imaris_Bridge.In", getArgument());
run("Unsharp Mask...");
call("Imaris_Bridge.Out", getArgument());
```

Sends processed volume back into Imaris Scene (over writing existing Volume object)



BITPLANE

an **Oxford Instruments** company



Imaris Open

Forum

File
Exchange

Imaris XT
Developer
Program

What Can the IO Platform Be Used for?

Organic FAQ

Meeting point
between Life and
Computer Scientists

A place for
collaborative work

Word-of-Mouth /
Reference Point

Centralized platform
for requesting and
sharing XTensions

Curated repository of
XTensions (citable,
unique and
permanent url)



Aaron
Ponti



Angela
Stathopoulos



Christopher
Wood



Jean-Yves
Tinevez



Jonas
Dorn



Josh
Thackray



Mario
Emmenlauer



Mary Cathleen
McKinney



Peter Beemiller



Lee Ling
(Sharon) Ong



Ricardo
Henriques



Richard
Alexander